

Marathon - Air walker

ART.NO 80200, 80201

Marathon is designed to increase stability and reduce edema. The orthosis has a built-in pump to achieve the desired compression. Shock-absorbing inner sole and non-slip outer sole with low profile promotes a natural walk, adjustable toe protection.

Indications: Indicated for stable foot and ankle fractures, severe sprains, soft tissue injuries, post-operative and post-cast removal. **Contraindications:** No known.

Material: 40% Glass fiber reinforced polyamide, 25% Polyester, 15% Polypropylene, 10% Thermoplastic rubber (TPR), 5% TPU, 5% EVA. Latex free.

Washing instructions: Liner should be hand washed only using warm water and detergent for delicate fabrics, repeatedly rinsed and line dried only. Do not use bleach or place in a dryer. Wipe plastic shell with damp cloth, you may use disinfectant or mild soap and wipe dry. Periodically check the condition of the Velcro to ensure they are free from grease and damage.

Accessories:

- Hygiene protection 10/box (Universal): 80219-00-00
- Socks, 2/box (Universal): 80215-30-00

Size	Shoe size	Art.no High	Art.no Ankle*
XS	<35	80200-00-30	80201-00-30
S	36-38	80200-00-40	80201-00-40
M	39-42	80200-00-50	80201-00-50
L	43-45	80200-00-60	80201-00-60
XL	>46	80200-00-70	80201-00-70

*Non stock item.



80200-High

Integrated pump in front to regulate the compression.



80201-Ankle*



MD Class I CE