

# StomiSoft - Stoma binder

ART.NO 50218, 50219

Straight elastic stoma binder for comfortable support of stoma and parastomal hernia. The binder has a non-elastic front. Easy to seal and adjust with functional handle and velcro fastening in the front. Velcros fastens on the entire surface of the non-elastic front material and makes the binder flexible in size and circumference. The binder can be adjusted to attain desired pressure.

Possible to cut hole in the non-elastic front. Holes can be cut in the desired size, adapted to the stoma bag. The front part does not fringe regardless of the shape and size of the hole, no sewing is required. Suitable to use with or without cut hole. If stonger compression is desired on the stoma bag, the binder closes so that the bag is underneath the non-elastic material. In order for the stoma bag to expand easily, the binder can be turned so that the bag is underneath the elastic material.

Reversible to suit stoma on both right and left side. Available in black, height 16 and 20 cm.

**Indications:** Stabilizing girdle for hernia/stoma/ostomy hernia. Prevents herniation. Fixation of stoma bandage. Mobilization after surgery.


**Contraindication:** incarcerated hernia.

**Material:** Polyester 49%, Polyamide 29%, Polyurethane 11%, Elastane 11%. Latex free.


**Washing instructions:** Machine wash at 40°C (106°F). Seal the Velcros. Use washing bag. Do not bleach. Do not dry clean. Do not iron. Do not tumble dry.

**Size/Measuring:** Measure circumference around intended placement. Pick the correct size according to the measurement table.


| Size | Circumference cm | Height 16 cm<br>Art.no | Height 20 cm<br>Art.no |
|------|------------------|------------------------|------------------------|
| XS   | 70-85            | 50219-00-30            | 50218-00-30            |
| S    | 80-95            | 50219-00-40            | 50218-00-40            |
| M    | 90-105           | 50219-00-50            | 50218-00-50            |
| L    | 100-115          | 50219-00-60            | 50218-00-60            |
| XL   | 110-125          | 50219-00-70            | 50218-00-70            |
| 2XL  | 120-135          | 50219-00-80            | 50218-00-80            |




Handle to easy put on and remove.



With hole cut for the stoma bandage.



With more pressure, non-elastic support over the stoma bandage.



With soft elastic support over the stoma bandage.



MD Class I CE



STANDARD  
100  
1976-298  
DTI